



# Clarius HD3 Vet Scanners

## Technical Specifications



### Clarius HD3 Vet Scanner Specifications

Model		Frequency	Max Depth	# Elements	Radius	Field of View	Pitch
<b>C3 HD3 Vet</b>	Convex	2–6 MHz	40 cm	192	45 mm	73°	300 µm
<b>C7 HD3 Vet</b>	Microconvex	3–10 MHz	18 cm	192	20 mm	112°	205 µm
<b>L7 HD3 Vet</b>	Linear	4–13 MHz	11 cm	192	N/A	38 mm	200 µm

### Imaging

#### Transmission

- 1 to 20 MHz waveforms
- Up to 20 continuous pulses
- Bi-polar output
- 10 to 70V peak-to-peak

#### Post-processing

- Adaptive speckle reduction
- Edge enhancement
- Persistence

#### Total Input Dynamic Range

- 160dB

#### Beamforming & Reception

- 8 parallel beamformers
- Synthetic aperture beamforming with virtual focal zones
- 60 MHz sampling rate @ 14 bits per channel

#### Automated Algorithms

- Time-gain-compensation (TGC)
- Frequency-depth adjustment
- Patient contact detection
- Needle enhancement
- Motion sensing
- Heart Rate

### Imaging Modes

<b>B-Mode</b>	Yes
<b>M-Mode</b>	Yes
<b>Power Doppler</b>	Yes
<b>Color Doppler</b>	Yes
<b>Pulsed-Wave Doppler</b>	Optional
<b>Elastography</b>	Optional
<b>Needle Enhance</b>	Optional

Available on L7 HD3 Vet

## Clinical Applications†

---

- C3 HD3 Vet**
- General Vet
  - General Cardiac
  - Abdomen
  - Cardiac
  - Bladder
  - MSK
  - Lung
  - Obstetrics

- 
- C7 HD3 Vet**
- General Vet
  - General Cardiac
  - Abdomen
  - Cardiac
  - Bladder
  - MSK
  - Lung
  - Obstetrics
  - Small Organs

- 
- L7 HD3 Vet**
- General Vet
  - MSK Equine
  - MSK
  - Lung
  - Nerve
  - Ocular
  - Vascular
  - Small Organs
- 

## Interface & Image Controls

---

- Depth
- Read zoom
- 3 TGC sliders or automated TGC
- Flip / mirror
- Freeze
- Color / power ROI
- Flow speed
- Doppler gate
- Doppler correction angle
- Doppler steer
- Baseline
- Invert

## Advanced Controls†

---

- Chroma
- Dynamic Range
- HD Zoom
- Trapezoidal
- Smoothing
- Penetration Mode

## Standard Configuration

---

- Scanner
- 1 Charger with global AC adapter

## Battery, Charging and Bootup

---

**Battery Life** ~60 min scanning

**Charge Time** ~90 min

**Bootup** Platform dependent,  
generally less than 30 sec

## Connectivity

---

**Wi-Fi** 802.11 a/b/g/n, dual band  
2.4GHz & 5GHz

**Bluetooth** Bluetooth low energy 4.1

## Warranty\*

---

**Included** 3 year limited warranty

**Optional** Clarius Care - 1/2/3 years

- Accidental damage
- Uptime
- RMA shipping
- Hospital theft
- Battery Service

† Some Presets available via Advanced Software Packages

\* [Click here](#) to see full terms and conditions

## Internally Optimized Parameters

Clarius internally optimizes the following parameters to ensure the scanner is easy to use:

<b>Frequency Range</b>	1 to 20 MHz
<b>Focal Zones Range</b>	1 to 10
<b>Compression Dynamic Range</b>	30 to 90 dB
<b>Reject</b>	Yes
<b>Sector Width Range</b>	50% to 100%
<b>Grey + Color Maps</b>	Yes
<b>Frame Rate</b>	Up to 30 FPS

## Mechanical

<b>Enclosure</b>	<ul style="list-style-type: none"><li>- Light weight magnesium alloy</li><li>- Durable</li><li>- IP67 rated for 1 meter immersion for 30 minutes</li></ul>
------------------	--

## SCANNER DIMENSIONS AND WEIGHT

<b>C3 HD3 Vet</b>	Dimensions: 146 x 76 x 32 mm Weight: 308 g
<b>C7 HD3 Vet</b>	Dimensions: 151 x 76 x 32 mm Weight: 289 g
<b>L7 HD3 Vet</b>	Dimensions: 147 x 76 x 32 mm Weight: 288 g

## CHARGER

<b>Input</b>	Wall power supply: 100-240 VAC, 50-60Hz Charger: 5 VDC, 3.2 A
<b>Output</b>	Wall power supply: 5 VDC, 3.2 A Charger: 5 VDC, 3.2 A

## Measurements and Calculations

### TOOLS

<b>Angle</b>	Yes
<b>Distance</b>	Yes
<b>Trace</b>	Yes
<b>Ellipse</b>	Yes
<b>Heart Rate</b>	Yes
<b>Time</b>	Yes
<b>Velocity</b>	Yes
<b>Volumes</b>	Yes → Manual/Automated

### CALCULATION PACKAGES

<b>Obstetrics</b>	HC, AC, CRL, GS, AFI, CxL, UA, FHR + Up to 12 Gestations
<b>Bladder</b>	Volume
<b>Abdomen</b>	Liver, Kidney, Spleen, Pancreas, GB
<b>Ocular</b>	ONSD, FB, CHLS
<b>Small Organs</b>	Thyroid, Testes
<b>Vascular</b>	Volume Flow
<b>Equine MSK</b>	SDFT, DDFT, ALDDFT, SL, ODSL, SDSL
<b>Cardiac</b>	LV EF (Simpsons), LV FAC, IVS, LVID, LVPW, TAPSE, MAPSE, EPSS, LVOT, HR

## Data Management

---

Local Export	JPG/PNG/DICOM/BMP
Cloud Export	Optional
DICOM Store	Optional
DICOM Worklist	Optional

## Security and Encryption

---

Wi-Fi Data Channel TLS 1.2

Bluetooth AES128 and RSA4096

## Standards Compliance

---

IEC 60601-1, Medical Electrical Equipment – Part 1: General requirements for basic safety and essential performance

IEC 60601-1-2, Medical Electrical Equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests

IEC 60601-2-37, Medical Electrical Equipment – Part 2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment

NEMA UD-2, Acoustic Output Measurement Standard For Diagnostic Ultrasound Equipment

NEMA UD-3, Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment

IEC 60601-1-12, Medical electrical equipment – Part 1-12: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems intended for use in the emergency medical services environment

## Cleaning

---

### Tested without adverse effects

- Accel® PREvention™ Wipes
- CaviWipes
- McKesson OPA/ 28 High-Level Disinfectant Solution
- MetriCide™ OPA Plus High-Level Disinfectant Solution
- Sani-Cloth® Plus Germicidal Disposable Cloth
- Tristel Trio Wipes System

FCC 47CFR Part 15, Radio frequency devices

ETSI EN 300 328 – Electromagnetic compatibility and Radio spectrum Matters (ERM)

ETSI EN 301 489-1 – Electromagnetic compatibility and Radio spectrum Matters (ERM)

ETSI EN 301 489-17 – Electromagnetic compatibility and Radio spectrum Matters (ERM)

ISO 10993-1, Biological evaluation of medical devices

IEC 60529, Degrees of protection provided by enclosures (IP Code)

IEC 62133, Secondary cells and batteries containing alkaline or other non-acid electrolytes – Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications

UN 38.3, Transport of dangerous goods – Classification procedures, test methods and criteria relating to class 9 – Lithium metal and lithium ion batteries

## About Us

Clarius Mobile Health was founded by experienced innovators who have played an instrumental role in the ultrasound industry. Our developers were the brains behind the first PC-based platform for ultrasound research. They also introduced the first touch screen ultrasound system with a simplified user interface.

We started with a simple mission: to enable more clinicians to use ultrasound to improve patient care. Thanks to the power of smart phones, advanced technology and decades of collective ultrasound experience, we've developed a high quality, Point-and-Shoot Ultrasound™ scanner that works with your smart device.

### Clarius Mobile Health

130-2985 Virtual Way  
Vancouver, BC V5M 4X7,  
Canada  
hello@clarius.com  
Phone: 1-778-800-9975