

## Transperineal Access System

*Revolutionising the methodology for obtaining prostate biopsies*

The PrecisionPoint™ Transperineal Access System is revolutionising prostate biopsies. The PrecisionPoint™ single use disposable device is a convenient, free hand innovation that enables transperineal biopsies without the need for expensive accessory equipment.


### The Benefits

- Enables urologists to easily transition to this new transperineal approach
- Decreases infections over the transrectal approach
- Allows systematic template-style sampling without the need for a template
- Avoids multiple perineal skin puncture
- Facilitates higher cancer detection rate
- Allows access to the anterior prostate
- Can be performed under local anesthesia



### Transrectal vs. Transperineal

By passing the biopsy needle through the perineum instead of the rectum, the significant risks associated with the transrectal approach are avoided.



# PrecisionPoint™

## Transperineal Access System

*Revolutionising the methodology for obtaining prostate biopsies*

The practitioner is not restricted to biopsy locations dictated by a grid configuration, since the technique is freehand. The practitioner can freely maneuver the ultrasound probe to align the access needle to target the desired locations with certainty for targeted and systematic biopsies.

It is designed to stabilize a Perineologic Access Needle on the anatomy as a means to properly guide the anatomical placement of biopsy needle into the prostate or other soft tissue through the use of the ultrasound-assisted guidance.

### About BXTAccelyon®

BXTAccelyon® is the sole distribution partner outside North America for the TGA approved PrecisionPoint™ Transperineal Access System.

### Contact

#### UK & Ireland:

**John Alden** john.alden@bxt-accelyon.com  
+44 (0) 7881 029 140

#### Benelux, Austria, Germany, Switzerland:

**Jos Groeskamp** jos.groeskamp@bxt-accelyon.com  
+31(0) 646808058

#### Australasia:

**Nicola Leavold** nicola.leavold@bxt-accelyon.com  
+61 409 458 852

**General enquiries:** customer@bxt-accelyon.com  
0800 014 8060  
+44 (0)1628 308 420

### Ordering Information

Product Ordering Code	Description
ACY AS 1031	PrecisionPoint™ Transperineal Access System - freehand biopsy device. 5 devices in a box PrecisionPoint™



facebook.com/bxtaccelyon



twitter.com/BXTAccelyon



linkedin.com/company/bxtaccelyon



youtube.com/channel/UChTU9L10CS5c1eWtGiMZ7gw