

## Gold Fiducial Markers- Now available to order with BXTAccelyon

- Electropolished for sharpness
- Siliconized for easy insertion
- Mirror polished I.D. for a smooth, easy seed deployment
- Echo-enhanced tipped for ultrasound visualization
- Centimetre markings for depth placement
- Rounded stylet for smooth gold marker delivery
- Foil marker to indicate bevel tip orientation



Gold Marker kits are available in 3 pre-loaded needles and 1 pre-loaded needle.

Gold Marker kits are supplied in square hub needle and are available in various marker sizes- 10 kits are supplied in each box.

<b>Gold Marker Kits of 3 pre-loaded needles</b>	
ACY GM 1033 1820	1 x 3mm Gold Marker- Gold Marker loaded in a pre-waxed 18G x 20cm needle
ACY GM 1233 1720	1.2 x 3mm Gold Marker- Gold Marker loaded in a pre-waxed 17G x 20cm needle
ACY GM 1053 1820	1 x 5mm Gold Marker- Gold Marker loaded in a pre-waxed 18G x 20cm needle
ACY GM 1253-1720	1.2 x 5mm Gold Marker- Gold Marker loaded in a pre-waxed 17G x 20cm needle

<b>Coupled Gold Marker Kits of 2 pre-loaded needles (four markers in total)</b>	
ACY GM 0852 SS20-1820	Gold Marker Kit, 18g, 20cm, 0.8 x 5mm, 2 needles, 2 Markers 20mm spacer with sleeve per needle

<b>Gold Marker Kits of 4 pre-loaded needles</b>	
ACY GM 1054-1820	1 x 5mm Gold Marker- 1.0mm Gold Marker loaded in a pre-waxed 18G x 20cm needle

**Other sizes and kits are available upon request.**

**For more information or to place your first Gold Fiducial Marker order please contact our customer service team on 0800 014 8060 or email: [customer@bxt-accelyon.com](mailto:customer@bxt-accelyon.com)  
[www.bxt-accelyon.com](http://www.bxt-accelyon.com)**