

## Clear Ultrasound Imaging You Can Trust

- ✓ No subscription required
- ✓ 3-year standard warranty
- ✓ Unlimited users
- ✓ Unlimited cloud storage



### Superior Image Quality

- High-quality images that are comparable to cart-based systems
- Imaging settings are optimized based on scanning application
- Optimal results from specialized scanners

### Reliable Portability

- Wireless freedom: no cables to trip over or break, and less risk of infections during procedures
- Stable connectivity from our proprietary multi-channel Wi-Fi architecture

### Easy to Use

- Point-and-Shoot Ultrasound™ technology removes the need to adjust knobs and switches
- Artificial intelligence adjusts parameters in real-time
- Compact and waterproof device that feels like standard scanners



Model	Frequency	Max Depth	Specialties
<b>C3 HD</b>	2–6 MHz	40 cm	EMED, EMS, OB/GYN, Primary Care, Urology
<b>L7 HD</b>	4–13 MHz	11 cm	Anesthesia, Breast, MSK, Pain, Surgery
<b>L15 HD</b>	5–15 MHz	7 cm	Anesthesia, Breast, MSK, Surgery
<b>L20 HD</b>	8–20 MHz	4 cm	Dermatology, MSK, Pediatrics, Plastic Surgery, Podiatry, Rheumatology, Vascular Access
<b>C7 HD</b>	3–10 MHz	18 cm	Neonatal, Pediatrics, Veterinary
<b>PA HD</b>	1–4 MHz	40 cm	Cardiac, Critical Care, EMED, EMS, OB/GYN
<b>EC7 HD</b>	3–10 MHz	15 cm	IVF / Reproductive, OB/GYN, Urology

#### Purchase Includes:

- M-Mode
- Power Doppler
- Color Doppler
- 1 hour remote in-service training
- Unlimited live telemedicine