

# Biopsy Suite<sup>®</sup>

## Transperineal Biopsy



Office-based urology chair with optional stepper-stabilizer system to perform transperineal prostate biopsy under local anesthesia

- ▶ Complete flexible system
- ▶ High quality
- ▶ Efficient and user-friendly

# Biopsy Suite<sup>®</sup> Accessories

## CHAIR

- Ref. code **400520** Urology chair, with electric height adjustment, pelvis section and back rest adjustment
- Ref. code **81717** Leg holders according to Goepel
- Ref. code **400066** Leg plate, for horizontal patient positioning, enlargement up to 1950 mm
- Ref. code **85223** Side rails L = 320 mm
- Ref. code **40023** Neck upholstery
- Ref. code **91617** Stainless steel disk
- Ref. code **40898** 4 double caster wheels Ø 75 mm
- Ref. code **40678** Potential equalisation socket
- Ref. code **40677** Socket at the base



## STAND / STEPPER

- Ref. code **UA-00-100-13B-0051**
- Articulate arm with six joints and tube connectors
  - Including universal OR-side rail adapter
  - Element diameter 50 mm
- Ref. code **BT-00-900-01-0194**
- Universal Brachytherapy-Stepper for all ER-probes
  - Prepared for two movement encoders for transperineal ultrasound probes
  - Please choose probe receptacle needed



## ACCESSORIES

- Ref. code **BT-00-900-02-0211**
- Single use, Disposable Template 18G
- Ref. code **BT-00-100-xx-0057**
- Reusable Template 17G



*Experts in diagnose and treatment technologies*

+31 (0)53 - 4836410  
info@jj-medtech.nl  
www.jj-medtech.nl

JJ-MT-BR-EN-001-A4-V2 Biopsy Suite